

Step-By-Step: Starting a C#.NET Data Acquisition solution using the Universal Library from scratch

Open C#.NET, and start a new Windows Application

Add a reference to MccDaq

- Go to your Solution Explorer window
- Expand the References folder (if it is not already expanded)
- Right mouse click on References folder, and select “Add Reference...”
- The Add Reference dialog box appears
- Scroll down to MccDaq, click on it, click on the Select button at the right
- Then click on the ‘OK’ button at the bottom of the dialog box

Add the namespace to the project

- Go to your Solution Explorer window
- click on Form1.vb then press F7 to activate the code window
- Locate the line “`using System.Data;`”
- Just after that add “`using MccDaq;`”

Add Project variables

- Locate the line “`publicclass Form1 : System.Windows.Forms.Form`”
- After the opening parenthesis add (paste in) the following code:

`//Here's where we declare our variables for the project`

```
private MccDaq.MccBoard DaqBoard;
```

```
private MccDaq.ErrorInfo ULStat;
```

```
private MccDaq.Range Range;
```

Step-By-Step: Starting a C#.NET Data Acquisition solution using the Universal Library from scratch

Add the Initialization to Form1 procedure

- Locate “`public Form1()`” procedure
- Locate “`InitializeComponent();`”
- Just after that add the following

```
// Initiate error handling

// activating error handling will trap errors like

// bad channel numbers and non-configured conditions.

// Parameters:

// MccDaq.ErrorReporting.PrintAll :all warnings and errors

// encountered will be printed

// MccDaq.ErrorHandling.StopAll :if an error is encountered,

// the program will stop
```

```
ULStat =
MccDaq.MccService.ErrHandling(MccDaq.ErrorReporting.PrintAll,
MccDaq.ErrorHandling.StopAll);
```

Add a new instance of a board object

- Go to your Solution Explorer window
- click on Form1.vb then press Shift + F7 to activate the Form Designer
- Double click on the form to activate the “`Form1_Load`” event
- Just after that add

Step-By-Step: Starting a C#.NET Data Acquisition solution using the Universal Library from scratch

```
int BoardNum;
```

```
BoardNum = 1; //board number assigned by InstaCal
```

```
DaqBoard = new MccDaq.MccBoard(BoardNum);
```

That's it. From here you can add code and procedures to things such as:

Please see the Universal Library function reference guide or C:\MCC\Unilib.chm for more information.

For example, you could add a command button and 2 labels to the form, then paste in the follow to the code window:

```
private void button1_Click(object sender, System.EventArgs e)
{
    float EngUnits;
```

```
System.UInt16 DataValue;
```

```
int Chan;
```

```
// Collect the data by calling AIn member function of MccBoard  
object
```

```
// Parameters:
```

```
// Chan :the input channel number
```

```
// Range :the Range for the board.
```

```
// DataValue :the name for the value collected
```

```
Range = Range.Bip10Volts; // select Bip10Volts (member of  
Range enumeration)
```

```
Chan = 0; // set input channel
```

```
ULStat = DaqBoard.AIn( Chan, Range, out DataValue);
```

```
if (ULStat.Value == MccDaq.ErrorInfo.ErrorCode.BadRange)
```

```
{
```

```
    MessageBox.Show( "Change the Range argument to one  
supported by this board.", "Unsupported Range",  
    MessageBoxButtons.OK);
```

```
    Application.Exit();
```

```
}
```

```
// Convert raw data to Volts by calling ToEngUnits (member  
function of MccBoard class)
```

```
ULStat = DaqBoard.ToEngUnits( Range, DataValue, out EngUnits);
```

```
label1.Text = DataValue.ToString();           // print the counts
```

```
label2.Text = EngUnits.ToString("F4") + " Volts"; // print the  
voltage
```

```
}
```

Note: This will only work for MCC hardware that has an A/D converter with a Bipolar 10 volt range.